National Center for Emerging and Zoonotic Infectious Diseases

Division of Foodborne, Waterborne and Environmental Diseases



Development of air-sampling and molecular detection methods for *Coccidioides* species

Anastasia Litvintseva

Mycotic Diseases Branch

Centers for Disease Control and Prevention

The Impact and Control of Valley Fever, NASEM Workshop November 17-18, 2022

Most people are infected by *Coccidioides* spp. from inhaling arthroconidia



No prevention aside from mitigating occupational risks

Utility of air-surveillance

- To better understand the geographic distribution
- To better understand the epidemiology:
 - Seasonality
 - Severe weather
 - Climate
 - Human factors
- To inform risk prediction models
- Generate "Cocci risk" forecasts
- Establish "Cocci warning" system



Technical challenges

- Identifying efficient air-sampling method
- Developing effective air-sampling strategy
- Developing sensitive and specific molecular detection method

Original Article

Molecular detection of airborne *Coccidioides* in Tucson, Arizona

Nancy A. Chow^{1,*}, Dale W. Griffin², Bridget M. Barker^{3,4,5}, Vladimir N. Loparev⁶ and Anastasia P. Litvintseva^{1,*}

¹Mycotic Diseases Branch, Centers for Disease Control and Prevention, Atlanta, Georgia, ²US Geological Survey, Coastal and Marine Science Center, St. Petersburg, Florida, ³Division of Pathogen Genomics, Translational Genomics Research Institute, Flagstaff, Arizona, ⁴Center for Microbial Genetics and Genomics, Northern Arizona University, Flagstaff, Arizona, ⁵Valley Fever Center for Excellence, University of Arizona, Tucson, Arizona and ⁶Biotechnology Core Facility Branch, Centers for Disease Control and Prevention, Atlanta, GA

Chow et al, 2017, Medical Mycology





Air samplers






High volume sampler (15 ft³/min)



Burkard mutli-stage liquid impinger (0.46 ft³/min)

Molecular detection

- Target unique transposable element in Coccidioides genome
- Approximately 70 copies per genome
- Not found in other fungal pathogens and related species
- Several variations of this assay, all target the same region



Air-samples were collected above known positive soils




Sampling ambient air



Sampling disturbed dust

Coccidioides DNA was detected in dust samples



Real-time PCR results from dust samples collected in the study

- Coccidioides DNA can be detected in disturbed dust samples using molecular methods
- Only high-volume samplers were successful
- Insufficient sampling time or methodology for detection in ambient air

Ambient Air Surveillance, 2016 Pilot in Arizona

- Daily collection of air-samples using high volume samplers
- Twenty-one sites around Phoenix, AZ
- Residual air filters shared from Sept 25-Nov 11, 2016



Air-filters collected in the pilot study








High-volume samplers




Air sampler

- Air-flow rate 100 L/min
- Collection time 24 h = 144,000 L of air sampled
- For comparison, a human inhales 8,000-9,000 L/day



Hydrophobic polytetrafluoroethylene (PTFE) filters

Modified DNA extraction










qPCR testing


STnested qPCR Limit of Detection

- Blank PTFE filters (without dust and other background) spiked with arthroconidia
- 10⁶ 1 arthroconidia per filter
- Able to detect one arthroconidium per filter
- Done in triplicates




Real-time PCR standard curve showing the correlation between added arthroconidia and Ct value

Daily prevalence of arthroconidia at 21 sites around Phoenix, AZ, September 25 - November 11, 2016



Temporal distribution of arthroconidia among 21 sites




https://academic.oup.com/mmy/article/58/4/552/5567250

Conclusions from the pilot

- Coccidioides can be successfully detected on PTFE air filters from ambient air
- Limit of detection: up to one spore on a filter in the laboratory
- Pilot testing results suggest uneven geographic prevalence and possible temporal signal

CocciWatch -Longitudinal air surveillance study









Same sites, 2017, 2018, 2019

Current datasets:

- Reported cases (ADHS)
- Air filters from Phoenix area (ADHS, DHS)
- Processed air filters (TGen, CDC)
- Weather data (TGen, NOAA)

Additional datasets we are working to include:

- Clinical testing data (positive/negative tests) (SQL)
- Valley fever cases in dogs
- Veterinarian testing data
















CocciWatch

- Filter collections from 23 sites (n = ~10K).
- Processing has focused on 11 of 23 sites, up to July 2019 (n = 5,055).





CocciWatch- Current Data




CocciWatch- Surveillance conclusions

- High variance in prevalence across space and time suggesting that risk is spatially and temporally variable.
 - Supported by ecology-Coccidioides must be introduced into the soil and then aerosolized.
- Local drivers influence site prevalence.
 - Specific weather patterns?
 - Soil disturbance?
 - Land cover around sites?








Summary

- Environmental surveillance is informative for understanding the epidemiology of Valley Fever
- "Patchy" distribution of arthroconidia in the air: Local conditions likely control the prevalence of arthroconidia in the air
- Air-monitoring can be informative for addressing specific questions and testing hypotheses; however, the results are site-specific

Possible applications

- To better understanding of the geographic distribution
 - To monitor emergence in new areas due to climate change
- To identify potential "hot spots" (for targeted vaccination)
- To identify human activities that release arthroconidia in the air to better understand and mitigate occupational risks

Acknowledgments







Lalitha Gade
Nancy Chow
Orion McCotter
Mitsuru Toda
Tom Chiller

Tanner PorterJolene Bowers
David Engelthaler

Dale Griffin





Thank you!

For more information, contact CDC 1-800-CDC-INFO (232-4636)
TTY: 1-888-232-6348 www.cdc.gov



The findings and conclusions in this report are those of the authors and do not necessarily represent the official position of the Centers for Disease Control and Prevention.

