

Strategies to Improve Characterization of Stem-Cell Based Cellular Products

National Academies of Sciences, Engineering, and
Medicine's

Forum on Regenerative Medicine

June 26, 2017


Steven R. Bauer, Ph.D.
Chief, Cellular and Tissue Therapies Branch
Division of Cellular and Gene Therapies
Office of Tissues and Advanced Therapies

Product Testing



Chemistry, Manufacturing, and Controls (CMC)

- Should ensure product safety
- Should ensure consistency of process and product
- Should (ideally) predict *in vivo* activity
- Is guided by detailed understanding of the manufacturing **process** and the **product**



Characterization


Multipotent Stromal Cells (MSCs) in Clinical Research



- **MSCs are diverse**
 - Characterization
 - Manufacturing
 - Source
- **Quality attributes?**
- **Relation to performance in clinical trials?**

MSC-Based Product Characterization for Clinical Trials: An FDA Perspective.
Cell Stem Cell:14, 2014 pg 141-145

MSCs and Product Characterization



?
=



**Manufactured
MSC**




Adherent in culture
CD73+, 90+, 105+
Lymphoid lineage negative




CBER/FDA MSC Consortium: Identification and correlation of MSC attributes with *in vivo* and *in vitro* assays of safety and efficacy


MSC Characterization




Puri Lab: genomics



McCright Lab: *in vivo*, *in vitro* models of wound repair




Bauer Lab: *in vitro* quantitative differentiation




PRODUCT CHARACTERISTICS

Moos Lab: gene expression, qRT-PCR, single cell PCR, NGS



CORRELATE CANDIDATE ATTRIBUTES WITH ASSAY OUTCOMES



Wei/Bauer Labs: *in vitro*, *in vivo* immunosuppression




Sung Lab: 3D *in vitro* tissue/organ microfluidic models




Hursh lab: epigenetics, karyotypes




Alterman Lab: proteomics





Consensus MSC Surface Markers Do Not Reveal Significant Differences



Adipogenic Potential Varies Between Cell Lines and Decreases with Passaging





Automated
microscopy-
based
quantification




Are Morphological Differences Associated With Different Biological Functions?

Lo Surdo JL et al.,
Cytotherapy, 2013.





MSC Morphology Predicts Osteogenic Activity

CULTURE-EXPANDED hMSCs




Marklein et al, Stem Cells 2016


Immunosuppressive Capacity Varies Between Cell Lines and with Culture Duration



↑ AUC ↓ Immunosuppressive Capacity



IFN γ -Stimulated MSC Morphology Predicts Immunosuppressive Activity



MSC Consortium Outcomes

- **Demonstrated that consensus markers do not predict functional biological heterogeneity of MSCs**
- **Identified cellular characteristics that predict relevant biological properties of MSC preparations**
 - Functionally relevant morphology
 - Gene expression
- **Published findings**
 - Sector Overview
 - Quantitation of Differentiation
 - Proteomics
 - Immunomodulation
 - Genomics
 - Genetic and Epigenetic Stability

Potential Applications (1)

- **Cell Source/Donor**
 - Screen samples for desired biological activity
- **Manufacturing**
 - Evaluate impact of manufacturing process
 - Tissue culture conditions and duration
- **Cell Characterization**
- **Identify Quality Attributes**
 - Activity/Potency
 - Quantitative Bioassays
 - Molecular markers correlated with bioassay outcomes
- **Guide cell enrichment techniques**

Potential Applications (2)

–Standards Development

- Quantitative bioassays
 - Osteogenesis
 - Adipogenesis
 - Immunosuppressive Capacity
 - Others?
 - »Angiogenesis
 - »Wound repair

MSC Consortium



Michail Alterman

- Samuel Mindaye
- Kristin Shultz-Kuszek
- Natalia Pripuzova



Steve Bauer

- Jessica Lo Surdo
- Heba Degheidy
- Ross Marklein
- Mike Mendicino
- Saniya Godil
- Matthew Klinker



Deb Hursh

- Patrick Lynch
- Yasmin Rovira-Gonzalez
- Brian Stultz
- Katie Steers



Kyung Sung

- Johnny Lam

Brent McCright

- Mandy Bush



Malcolm Moos

- Sema Rosinbum
- Jennifer Mateshaytis
- Elaine Thompson
- Alexa Bianchi



Raj Puri

- Ian Bellayr
- Jennifer Catalano



Cheng-Hong Wei

- Cristina Nazarov
- Matthew Klinker



- CBER/OTAT/DCGT
- FDA Targeted Research Funds
- BARDA

Contact Information

- **Steven R. Bauer, Ph.D.**

steven.bauer@fda.hhs.gov

- **Regulatory Questions:**

OTAT Main Line – 240 402 8190

Email: OTATRPMS@fda.hhs.gov and

Lori.Tull@fda.hhs.gov



*FDA Headquarters
Federal Research Center at White Oak
10903 New Hampshire Avenue
Silver Spring, MD 20993-0002*

- **OCTGT Learn Webinar Series:**

<http://www.fda.gov/BiologicsBloodVaccines/NewsEvents/ucm232821.htm>

- **CBER website:** www.fda.gov/BiologicsBloodVaccines/default.htm

- **Phone:** 1-800-835-4709 or 240-402-8010

- **Consumer Affairs Branch:** ocod@fda.hhs.gov

- **Manufacturers Assistance and Technical Training Branch:** industry.biologics@fda.gov

- **Follow us on Twitter:** <https://www.twitter.com/fdacber>

