



ReadiMask[®]

*Protect yourself
Protect others*

PPE for the Public
January 26th 20216



Presentation

- ✓ A demonstration
- ✓ ReadMask History
- ✓ Respirators should...
- ✓ ReadMask design and adaptation
- ✓ Q&A



A demonstration of our strapless, NIOSH approved N95 respirator

Fits All face shapes and sizes

Guaranteed to pass QNFT & QLFT

Source Control

Comfort, Comfort, Comfort



“We love, love, love these masks!”

*Protect yourself
Protect others*



The ReadiMask History (the short version)

1998 – WTC

2010 – Hello NIOSH ☺


2019 – N95 Approval

2020 – SARS-CoV-2

Avery

DLA

AONS, A22Y



Respirators should be comfortable and tolerable

A man does to the doctor...



“Clear vision is overrated.”




“I never touch my mask!”



*Protect yourself
Protect others*



Respirators should seal



*Protect yourself
Protect others*



Respirators should be age and body size appropriate



Respirators should be trusted



*Protect yourself
Protect others*



Respirators should be tested on people as part of certification



FIT TEST REPORT
12/4/2020

ID NUMBER	41365	CUSTOM1	
LAST NAME	FOSTER	CUSTOM2	
FIRST NAME	JUDY	CUSTOM3	
COMPANY	RIT	CUSTOM4	
LOCATION	EH9		
TEST DATE	12/4/2020 12:19	PORTACOUNT S/N	8038103812
DUE DATE	12/4/2021	N95 COMPANION	Y
RESPIRATOR	READIMASK N1902S N95 [169]	PROTOCOL	OSHA 29CFR1910.134
MANUFACTURER	READIMASK	PASS LEVEL	100
MODEL	N1902S	APPROVAL	
MASK STYLE	N95	EFFICIENCY-95%	True
MASK SIZE	small		

EXERCISE	DURATION (sec.)	FIT FACTOR	PASS
NORMAL BREATHING	60	200+	Y
DEEP BREATHING	60	200+	Y
HEAD SIDE TO SIDE	60	200+	Y
HEAD UP AND DOWN	60	200+	Y
TALKING	60	200+	Y
GRIMACE	15	Excl.	
BENDING OVER	60	200+	Y
NORMAL BREATHING	60	200+	Y
OVERALL FF		200+	Y

FIT TEST OPERATOR _____ DATE _____
JUDY FOSTER

NAME _____ DATE _____
JUDY FOSTER

Note: WEARS GLASSES

Respirator Fit Test Card
Name: JUDY FOSTER Test Date: 12/4/2020
ID: 41365 Next Test Date: 12/4/2021

Respirator	Results
Mfg: READIMASK	Overall FF: 200+
Model: N1902S	FF Pass Level: 100
Style: N95	Pass: Y
Size: small	Operator: JUDY FO...

Protocol: OSHA 29CFR1910.134
Fit Test Method: CNFT using TSI PortaCount
*** Your company contact information here ***

ReadiMask: Guaranteed to fit

*Protect yourself
Protect others*




How it's made

Universal Fit, Easy to customize

Layers:

1. Aperture film
2. Meltblown
3. Spunbond
4. Adhesives:
 - Silicone
 - Acrylic



Protect yourself
Protect others



Where it's made

HQ – New Jersey

MFG – MD & OH

Materials – DE, OH, MN

Berry compliant



RadiMask with integrated eye shield



Samples and Meeting Request

Email – Jschwind@ReadiMask.com

Website – ReadiMask.com

LinkedIn

*Protect yourself
Protect others*



Q & A



*Protect yourself
Protect others*

John Schwind
Jschwind@ReadiMask.com
732-616-6781

